


Clearview **Simplify D-dimer**

Fingerprick blood
 Kapillarblut aus der Fingerbeere
 Fingerprikning
 Colecta de Sangre por Punción Digital
 Sang du bout du doigt
 Fingerstick Blood
 Sangue dal polpastrello
 Vingerprikbloed
 Blod fra fingertupp
 Sangue por picada do dedo
 Kapillärblod med fingerstick

Venous whole blood
 Venöses Vollblut
 Venøst helblod
 Sangre Venosa Total
 Sang veineux total
 Φλεβικό ολόκληρο αίμα
 Sangue intero venoso
 Veneus vol bloed
 Venøst fullblod
 Sangue Venoso Total
 Venøst helblod

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Analyse	Sensitivitet	Specificitet	NPV
Clearview Simplify D-dimer	100% (90.0-100%)*	52.9% (41.8-63.8%)*	100% (92.1-100%)*
Vidas® D-dimer New	100% (90.7-100%)*	48.8% (37.6-60.1%)*	100% (91.2-100%)*

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Analyse	Sensitivitet	Specificitet	NPV
Clearview Simplify D-dimer	100% (92.5-100%)*	47.9% (43.3-52.6%)*	100% (98.4-100%)*
Vidas® D-dimer New	100% (92.5-100%)*	48.8% (44.1-53.4%)*	100% (98.4-100%)*
STA® Liatest D-DI	100% (92.5-100%)*	47.5% (42.7-52.3%)*	100% (98.2-100%)*

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www.clearview.com

product.support@invmed.com

© 2009 Inverness Medical. All rights reserved. Clearview is a trademark of Inverness Medical group of companies.

British Biocell International Limited, Dundee, DD2 1NH, UK

Bibliografia/Referenties/Referencias/References/Henvisninger/Referenser/Literatur/Références/Referanser/Referências/Παρατιμπές

1. Gaffney, P.J. D-dimer History of Discovery, Characterisation and Utility of this and other Fibrin Fragments. *Fibrinolysis* 7 Suppl 2:2-8; 1993.
2. Lane, D.A. et al. Characterisation of Serum Fibrinogen and Fibrin Fragments Produced During Disseminated Intravascular Coagulation. *Haematology*. 40: 609-615; 1978.
3. Rylatt, D.B. et al. An Immunoassay for Human D-dimer using Monoclonal Antibodies. *Thromb. Res.* 31: 767-778; 1983.
4. CLSI Publication H18-A3 - Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition; 2004.
5. Wells, P.S. Evaluation of D-dimer in the Diagnosis of Suspected Deep-Vein Thrombosis. *N. Engl. J. Med.* 349 (13): 1227-1235; 2003.
6. Keeling, D.M. et al. The Haemostasis and Thrombosis Task Force of the British Committee for Standards in Haematology. The diagnosis of deep vein thrombosis in symptomatic outpatients and the potential for clinical assessment and D-dimer assays to reduce the need for diagnostic imaging. *Br. J. Haematol.* 124(1): 15-25; 2004.
7. Bick, R.L. et al. Diagnostic Efficacy of the D-dimer assay in Disseminated Intravascular Coagulation (DIC) *Thromb. Res.* 65:785-790; 1992.
8. Bick, R.L. et al. Disseminated Intravascular Coagulation: Objective Clinical and Laboratory Diagnosis, Treatment, and Assessment of Therapeutic Response. *Semin. Thromb. Hemost.* 22(1): 69-88; 1996.
9. Scarvelis, D and Wells, P.S. Diagnosis and Treatment of Deep Vein Thrombosis. *Can. Med. Assoc. J.* 175 (9):1087-92; 2006.
10. Subramanian, R.M. et. al. Does an Immunochromatographic D-dimer exclude acute lower limb deep venous thrombosis? *Emer. Med. Austral.* 18: 457-463; 2006.
11. Runyon, M.S. et. al. Comparison of the Simplify D-dimer assay performed at the bedside with a laboratory based quantitative D-dimer assay for the diagnosis of pulmonary embolism in a low prevalence emergency department population. *Emerg. Med. J.* 25:70-75; 2008.
12. Ginsburg, J.S. et. al. Sensitivity and specificity of a rapid whole-blood assay for D-dimer in the diagnosis of pulmonary embolism. *Ann. Intern. Med.* 129(12), 1006-11; 1998.
13. Hunt, F.A. et al. Serum Cross-Linked Fibrin (XDP) and Fibrinogen/Fibrin Degradation Products (FDP) in Disorders Associated with Activation of the Coagulation or Fibrinolytic Systems. *Br. J. Haematol.* 60: 715-722; 1985.
14. Smith, R.T. et al. Fibrin Degradation Products in the Post-Operative Period- Evaluation of a New Latex Agglutination Method. *AJCP.* 60: 644-647; 1973.
15. Nolan, T.E. et al. Maternal Plasma D-dimer Levels in Normal and Complicated Pregnancies. *Obstetrics & Gynecology.* 81(2): 235-238, 1993.
16. Cini M, Legnani C, et al. A New Rapid Beside Assay for D-dimer Measurement (Simplify D-dimer) in the Diagnostic Work-up for Deep Vein Thrombosis. *JTH Vol 1* (12): p 2681; 2003.
17. Toulon P, Meyniard O. D-dimer testing in patients with suspected pulmonary embolism. Performance of a new rapid qualitative solid-phase immuno-chromatography assay (Simplify D-dimer) and comparison with two automated quantitative assays. *Blood.* Vol 102, 11 (1b-514b) Abstract #4149; 2003.

Solo per uso diagnostico in vitro/Uitsluitend voor In Vitro Diagnostiek/Exclusivamente para diagnóstico in vitro/For In Vitro Diagnostic Use/Til in vitro-diagnose brug/För diagnostisk användning in vitro/Für die In-Vitro-Diagnostik/Destiné à un usage diagnostique in vitro/Til bruk ved in vitro diagnostikk/Para uso em Diagnóstico In Vitro/**Μόνο για διαγνωστική χρήση in vitro**

Numero di partita/Lotnummer/Número de Lote/Lot Number/Lot nummer/Satsnummer/Chargennummer/Numéro de lot/Batch nummer/Número de lote/**Αριθμός παρτίδας**

Tampone/Buffer/Tampón/Buffer/Buffer/Buffer/Puffer/Tampon/Buffer/Soluto-tampão/**Ρυθμιστικό διάλυμα**

Consultare le istruzioni per l'uso/De gebruiksaanwijzing raadplegen/Consulte las instrucciones para el Uso/Consult Instructions for Use/Se brugsanvisning/Se bruksanvisningen/Gebrauchsanleitung konsultieren/Lire le mode d'emploi/Les bruksanvisningen/Consultar as Instruções para utilização/**Βλέπε τις Οδηγίες χρήσης**

Data di scadenza/Te gebruiken voor/Utilizar Antes de/Fecha de Caducidad/Use By/Expiry Date/Udløbsdato/Usar até/Använd före/Utgångsdatum/Verwendbar bis/A utiliser avant/Date de péremption/Anv. inden/Brukes ved/Prazo de Validade/**Χρήση ως/ Ημερομηνία λήξης**

Conservare a temperature comprese tra 2°C e 25°C/Opstaan bij 2-25°C/Almacenar a 2-25°C/Store at 2-25°C/Opbevares ved 2-25°C/Förvaras vid 2-25°C/Bei 2-25°C lagern/A conserver à une température comprise entre 2 et 25°C/Lagres ved 2-25°C/Conservar a 2°C-25°C/**Φυλάξτε στους 2-25°C**

Dispositivo per test/Testtoestel/Dispositivo para Prueba/Test Device/Testenhet/Analyseenhet/Testvorrichtung/Dispositif d'examen/Testapparat/Dispositivo para testes/**Δοκιμαστική συσκευή**

Produttore/Fabrikant/Fabricante/Manufacturer/Producent/Tillverkare/Hersteller/Fabricant/Produsent/Fabricante/**Παρασκευαστής**

Attenzione, vedere le istruzioni per l'uso/Voorzichtig, raadpleeg bijgevoegde documenten/Atención, ver instrucciones de uso/Caution, consult accompanying documents/Forsigtig se brugsanvisning/Försiktighet, se handhavandebeskrivningen/Achtung, Begleitdokumente beachten/Attention voir notice d'instructions/Advarsel: Les de vedlagte dokumentene/Atenção, consulte a documentação incluída/**Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα /**

Numero di catalogo / Catalogus nummer / Número de Catálogo / Catalogue Number / Katalognummer / Katalognummer / Katalognummer / Numéro de référence dans le catalogue / Katalognummer / Número de Catálogo / **Αριθμός καταλόγου**