
Clearview **Simplify D-dimer**

Fingerprick blood

Kapillarblut aus der Fingerbeere

Fingerprikning

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www.clearview.com

product.support@invmed.com

© 2009 Inverness Medical. All rights reserved. Clearview is a trademark of Inverness Medical group of companies.

British Biocell International Limited, Dundee, DD2 1NH, UK

Bibliografia/Referenties/Referencias/References/Henvisninger/Referenser/Literatur/Références/Referanser/Referências/Παρατιμπές

1. Gaffney, P.J. D-dimer History of Discovery, Characterisation and Utility of this and other Fibrin Fragments. *Fibrinolysis* 7 Suppl 2:2-8; 1993.
2. Lane, D.A. et al. Characterisation of Serum Fibrinogen and Fibrin Fragments Produced During Disseminated Intravascular Coagulation. *Haematology*. 40: 609-615; 1978.
3. Rylatt, D.B. et al. An Immunoassay for Human D-dimer using Monoclonal Antibodies. *Thromb. Res.* 31: 767-778; 1983.
4. CLSI Publication H18-A3 - Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition; 2004.
5. Wells, P.S. Evaluation of D-dimer in the Diagnosis of Suspected Deep-Vein Thrombosis. *N. Engl. J. Med.* 349 (13): 1227-1235; 2003.
6. Keeling, D.M. et al. The Haemostasis and Thrombosis Task Force of the British Committee for Standards in Haematology. The diagnosis of deep vein thrombosis in symptomatic outpatients and the potential for clinical assessment and D-dimer assays to reduce the need for diagnostic imaging. *Br. J. Haematol.* 124(1): 15-25; 2004.
7. Bick, R.L. et al. Diagnostic Efficacy of the D-dimer assay in Disseminated Intravascular Coagulation (DIC) *Thromb. Res.* 65:785-790; 1992.
8. Bick, R.L. et al. Disseminated Intravascular Coagulation: Objective Clinical and Laboratory Diagnosis, Treatment, and Assessment of Therapeutic Response. *Semin. Thromb. Hemost.* 22(1): 69-88; 1996.
9. Scarvelis, D and Wells, P.S. Diagnosis and Treatment of Deep Vein Thrombosis. *Can. Med. Assoc. J.* 175 (9):1087-92; 2006.
10. Subramanian, R.M. et. al. Does an Immunochromatographic D-dimer exclude acute lower limb deep venous thrombosis? *Emer. Med. Austral.* 18: 457-463; 2006.
11. Runyon, M.S. et. al. Comparison of the Simplify D-dimer assay performed at the bedside with a laboratory based quantitative D-dimer assay for the diagnosis of pulmonary embolism in a low prevalence emergency department population. *Emerg. Med. J.* 25:70-75; 2008.
12. Ginsburg, J.S. et. al. Sensitivity and specificity of a rapid whole-blood assay for D-dimer in the diagnosis of pulmonary embolism. *Ann. Intern. Med.* 129(12), 1006-11; 1998.
13. Hunt, F.A. et al. Serum Cross-Linked Fibrin (XDP) and Fibrinogen/Fibrin Degradation Products (FDP) in Disorders Associated with Activation of the Coagulation or Fibrinolytic Systems. *Br. J. Haematol.* 60: 715-722; 1985.
14. Smith, R.T. et al. Fibrin Degradation Products in the Post-Operative Period-Evaluation of a New Latex Agglutination Method. *AJCP.* 60: 644-647; 1973.
15. Nolan, T.E. et al. Maternal Plasma D-dimer Levels in Normal and Complicated Pregnancies. *Obstetrics & Gynecology.* 81(2): 235-238, 1993.
16. Cini M, Legnani C, et al. A New Rapid Beside Assay for D-dimer Measurement (Simplify D-dimer) in the Diagnostic Work-up for Deep Vein Thrombosis. *JTH Vol 1* (12): p 2681; 2003.
17. Toulon P, Meyniard O. D-dimer testing in patients with suspected pulmonary embolism. Performance of a new rapid qualitative solid-phase immuno-chromatography assay (Simplify D-dimer) and comparison with two automated quantitative assays. *Blood.* Vol 102, 11 (1b-514b) Abstract #4149; 2003.

Solo per uso diagnostico in vitro/Uitsluitend voor In Vitro Diagnostiek/Exclusivamente para diagnóstico in vitro/For In Vitro Diagnostic Use/Til in vitro-diagnose brug/För diagnostisk användning in vitro/Für die In-Vitro-Diagnostik/Destiné à un usage diagnostique in vitro/Til bruk ved in vitro diagnostikk/Para uso em Diagnóstico In Vitro/**Μόνο για διαγνωστική χρήση in vitro**

Numero di partita/Lotnummer/Número de Lote/Lot Number/Lot nummer/Satsnummer/Chargennummer/Numéro de lot/Batch nummer/Número de lote/**Αριθμός παρτίδας**

Tampone/Buffer/Tampón/Buffer/Buffer/Buffer/Puffer/Tampon/Buffer/Soluto-tampão/**Ρυθμιστικό διάλυμα**

Consultare le istruzioni per l'uso/De gebruiksaanwijzing raadplegen/Consulte las instrucciones para el Uso/Consult Instructions for Use/Se brugsanvisning/Se bruksanvisningen/Gebrauchsanleitung konsultieren/Lire le mode d'emploi/Les bruksanvisningen/Consultar as Instruções para utilização/**Βλέπε τις Οδηγίες χρήσης**

Data di scadenza/Te gebruiken voor/Utilizar Antes de/Fecha de Caducidad/Use By/Expiry Date/Udløbsdato/Usar até/Använd före/Utgångsdatum/Verwendbar bis/A utiliser avant/Date de péremption/Anv. inden/Brukes ved/Prazo de Validade/**Χρήση ως/ Ημερομηνία λήξης**

Conservare a temperature comprese tra 2°C e 25°C/Opstaan bij 2-25°C/Almacenar a 2-25°C/Store at 2-25°C/Opbevares ved 2-25°C/Förvaras vid 2-25°C/Bei 2-25°C lagern/A conserver à une température comprise entre 2 et 25°C/Lagres ved 2-25°C/Conservar a 2°C-25°C/**Φυλάξτε στους 2-25°C**

Dispositivo per test/Testtoestel/Dispositivo para Prueba/Test Device/Testenhed/Analyseenhet/Testvorrichtung/Dispositif d'examen/Testapparat/Dispositivo para testes/**Δοκιμαστική συσκευή**

Produttore/Fabrikant/Fabricante/Manufacturer/Producent/Tillverkare/Hersteller/Fabricant/Produsent/Fabricante/**Παρασκευαστής**

Attenzione, vedere le istruzioni per l'uso/Voorzichtig, raadpleeg bijgevoegde documenten/Atención, ver instrucciones de uso/Caution, consult accompanying documents/Forsigtig se brugsanvisning/Försiktighet, se handhavandebeskrivningen/Achtung, Begleitdokumente beachten/Attention voir notice d'instructions/Advarsel: Les de vedlagte dokumentene/Atenção, consulte a documentação incluída/**Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα /**

Numero di catalogo / Catalogus nummer / Número de Catálogo / Catalogue Number / Katalognummer / Katalognummer / Katalognummer / Numéro de référence dans le catalogue / Katalognummer / Número de Catálogo / **Αριθμός καταλόγου**